



ACCESSORIES

The Expanded Duty Cotton Mouth™

Impression Trays

Cheek Retractors

Lip Retractors

Facebows

Photo Mirrors

Cephalometric Accessories

Curing Lights

Patient Care Kits

Plastic zip pouch with essential oral care products.

Patient Care Kits - Blue (10 per pk)

#700-207

- Disposable Toothbrush
- Travel Toothbrush
- Dental Floss
- Sand Timer

Patient Care Kits - Pink (10 per pk)

#700-208

- Disposable Mirror
- Interdental Brush
- Floss Threaders
- Patient Relief Wax



The Expanded Duty Cotton Mouth™

The Expanded Duty Cotton Mouth™ provides a dry field for numerous dental and orthodontic procedures, while also affording your patient total comfort. The molar pad extenders allow the bite opening to be increased for maximum access.

Pad Extender-5mm (4/pk)

#900-204

Aspirator Tip-long (2/pk)

#900-203

Cotton Mouth Kit

- (1) Bite Block
 - (4) 5mm Extender Pad
 - (4) 7mm Extender Pad
 - (2) Aspirator Tip-short
 - (2) Aspirator Tip-long
 - (25) Disposable Tubes
- #900-200**

Pad Extender-7mm (4/pk)

#900-205

Bite Block (1/pk)

#900-201

Aspirator Tip-short (2/pk)

#900-202

Disposable Tubes (25/pk)

#900-206



Arch Markers

Thick wax tips upon sturdy, wooden dowels to provide precise marking, high visibility and superior adhesion. Disposable to prevent cross contamination.



100 Per Bag

#800-446



Brace Guard Squares

The ultimate in bracket protection, simply pinch over the offending bracket, pushing material under and behind the bracket, forming a gripping configuration (mechanical retention). Special food grade silicone wax will not melt or harden.

80 Per Box

#700-204



Patient Relief Wax

Ideal for fast, temporary relief from fixed brackets and appliances.

50 Per Box

#700-104



Access Cheek Retractor

Unique bow design allows adjustment to ensure comfortable fit. Offers excellent accessibility to bicuspid and molar area. May be autoclaved or dry heat sterilized.

Adult (1/pk)
#600-752

Pedo (1/pk)
#600-751



Expando Cheek Retractor

Allows full view of posterior teeth and is self-retaining. Excellent for photography. Can be cold sterilized or autoclaved up to 250°F.

Assortment (1 of ea.)
#600-049

Adult (2/pk)
#600-050

Teen (2/pk)
#600-052

Pedo (2/pk)
#600-051



Photo Lip Retractor -45°

Cold sterilization is recommended for extended life, but can be autoclaved up to 250°F.

Large (1/pk)
#600-034

Small (1/pk)
#600-035



Lip Retractor

Hand held retractor for upper and lower arch photography. Cold sterilization is recommended for extended life, but can be autoclaved up to 250°F.

#600-753



Retainer Case

Retainer Cases Mixed 24/pk - Sparkle
#700-028

Retainer Cases Mixed 24/pk - Pearl
#700-029



Multi-Adjustable Protraction Mask

An ideal ergonomic design to obtain maximum class III compliance. It is fully adjustable to each patient's profile and features a lightweight, durable design for maximum patient comfort.

Replacement Pads
#550-104

Blue (1/pk)
#550-100

Red (1/pk)
#550-101

Pearl (1/pk)
#550-102

Clear (1/pk)
#550-103



*Refer to page 56 for Intraoral Elastics.

Lip Bumper (Step down leg)

2/pk	Measurement	SKU
Size 2	86mm	#271-342
Size 3	92mm	#271-343
Size 4	98mm	#271-344
Kit	Size 2-4 (1 each)	#271-340



Side View



Cushioned Lip Bumper


2/pk	Measurement	SKU
Size 1	80mm	#271-351
Size 2	86mm	#271-352
Size 3	92mm	#271-353
Size 4	98mm	#271-354
Size 5	104mm	#271-355
Size 6	110mm	#271-356
Size 7	116mm	#271-357
Kit	Size 1-7 (1 each)	#271-350



Stainless Steel Photo Mirrors

These Intraoral Stainless Steel Photo Mirrors are highly polished, double-ended mirrors that provide the perfect imaging solution for every patient. The highly polished surfaces of Ortho Arch's stainless steel mirrors ensure high reflectivity. All of our intraoral stainless steel double ended mirrors are set at a 30° angle.


Adult Occlusal/Child Occlusal
#801-020




Buccal Lateral
#801-021



Child Occlusal/Child Buccal
#801-026




Adult Occlusal/Child Occlusal
#801-022




Adult Buccal/Adult Occlusal
#801-023



Lingual Upper/Lingual Lower
#801-024



Child Buccal/Adult Buccal
#801-025



Millimeter Ruler

10/pk




#700-064

Cephalometric Ruler/Protractor

A combination tooth template and protractor for accurate cephalometric analysis.

#700-069



Cephalometric Tracing Paper

Smooth, matte surface and transparent quality gives extra clean, sharp lines when tracing Radiographs.

8" x 10" sheets
100 Sheets per pad
#700-058



Disposable Impression Trays

Molded of high density plastic, these impression trays are excellent for mechanical retention of alginates without need for adhesives. The six color-coded sizes ensure an accurate fit.

1	Upper (10/pk) #700-220	2	Upper (10/pk) #700-224	3	Upper (10/pk) #700-228
	Lower (10/pk) #700-222		Lower (10/pk) #700-226		Lower (10/pk) #700-230
4	Upper (10/pk) #700-232	5	Upper (10/pk) #700-236	6	Upper (10/pk) #700-240
	Lower (10/pk) #700-234		Lower (10/pk) #700-238		Lower (10/pk) #700-242
	Assortment Kit (36/pk) #700-250				



Sport Mouthguard

Specially crafted with the thought of athletes with braces. Protects both the maxillary and mandibular teeth. Protects from lacerations to the lips and cheeks. The perfect mouthguard suited for all sports.

- #600-780 Sport Mouthguard - Clear (10)
- #600-781 Sport Mouthguard - Blue (10)
- #600-782 Sport Mouthguard - Silver (10)
- #600-783 Sport Mouthguard - Green (10)
- #600-784 Sport Mouthguard - Yellow (10)



Bracketed Typodont

1/pk	Size	SKU
Nustar Gen 2	7x7	#700-357NU
Horizon Mini	7x7	#700-357HZ
Micro Twin	5x5	#700-357MC
Nova C Ceramic	5x5	#700-357NC
Nova S Sapphire	5x5	#700-357NS
AccuTwin	5x5	#700-357AT
Conlin	5x5	#700-357CX
Clear-Lock	5x5	#700-357CL
Solo Twin	5x5	#700-357ST
D2 Appliance	-	#700-357D2
Class I Dentition	7x7	#700-357



ACCESSORIES

Typodont Class I Dentition (7x7)

#700-357



Clear Articulated Typodont

Class I Dentition

This clear typodont has removable teeth with anatomical roots. The clear flexible feature is excellent for patient education. The typodont is ideal for bonding brackets and placing molar bands.

Natural color teeth

#700-353

White teeth

#700-353W



#700-353

Model with Removable Teeth

Class I Dentition

An economical option for patient education or bonding brackets.

#700-340



Teeth Models

Anatomical teeth (second molar to second molar). Brackets may be bonded to this set of teeth.

#700-354



Super Fine Diamond Burs

These Super Fine Diamond Burs are ideal for interproximal tooth reduction, and contouring. Comes in a regular and long handle shank. Available in regular and long shank friction grip.

Short (5/pk)

#700-466

Long (5/pk)

#700-467



Diamond Disc

Used for contouring and interproximal reduction. Ideal for creating intra-arch space or opening contact points. Use with low speed handpieces only. They are heat sterilizable up to 340F/170C or autoclavable up to 250F/121C.

Double Sided- mandrel: 22mm, disc: 0.17mm (2/pk with mandrel)

#700-461

Single Sided- mandrel: 22mm, disc: 0.10mm (2/pk with mandrel)

#700-460

Double Sided See Thru- mandrel: 22mm, disc: 0.20mm (2/pk with mandrel)

#700-462



Honeycomb Diamond Disk

The Honeycomb design allows for smooth rotation, improved visibility, greater flexibility, better contouring, and less gumming of the disc.

#700-481



Single Sided Diamond Strips

.013 Single Sided Diamond Strips - Blue (10/pk) Medium

#700-482

.010 Single Sided Diamond Strips - Red (10/pk) Fine

#700-483

.08 Single Sided Diamond Strips - Yellow (10/pk) Ultra Fine

#700-484



Diamond Disc Safety Guard

The Diamond Disc Safety Guard can be used with 22mm E Style hand pieces and is autoclavable up to 275F/135C.

#700-463



Interproximal Reduction Gauge

Measures .10 mm, .20mm, .25mm, .30mm, .40mm and .50mm.

#800-139



Silicone Color Coded Instrument Rings

Assorted Colors.

80 per pack.

#700-050



REDI SET LED Curing Light

The LED Curing Light has a variable time setting from 5 to 40 seconds, operable up to 200 applications per charge. The optical wave length is from 420 - 480nm achieves a blue light intensity of 1600 - 1800mW/cm² Includes a 1 year warranty.

#800-511

- Delivers dependable curing for both stainless steel and ceramic brackets
- With a full charge can provide up to 200 applications
- Three work modes: Full power, ramping-up and pulse
- Variable settings from 5 - 40 seconds
- 420nm - 480nm wave length
- Light intensity of 1,600 mW/cm² to 1,800 mW/cm²
- Carries detection light as an added feature
- Includes a 1 year manufacturer warranty



Mini Curing Light

The Mini Curing Light has a spectrum range of 450-490nm, high light intensity, is autoclavable, contains a 1,000 hour 'solid state' LED bulb, and a life time warranty. Requires AAAA batteries.



#800-510

QWIKSTRIP

Finally, an accurate IPR system at your fingertips.
 The Qwikstrip abrasive strip system provides the ultimate IPR solution.
 Simply follow the color-coded sequence guide to achieve precise interproximal reduction.



- #700-400 QWIKSTRIP MASTER IPR KIT - (90: 5EA OF ALL GRITS & 10 SERRATED)
 - #700-401 QWIKSTRIP STARTER IPR KIT - SINGLE SIDED (10: 2 EA GRIT)*
 - #700-402 QWIKSTRIP STARTER IPR KIT - CURVED (10: 2 EA GRIT)*
 - #700-403 QWIKSTRIP STARTER IPR KIT - DOUBLE SIDED (10: 2 EA GRIT)*
 - #700-404 QWIKSTRIP STARTER IPR KIT - ASSORTED (10: 1 EA GRIT IN SGL & DBL-SIDED)*
- * Does NOT include Serrated or XX Coarse Strips.



- #700-411 QWIKSTRIP IPR YELLOW 0.07MM - SINGLE-SIDED (10)
- #700-421 QWIKSTRIP IPR YELLOW 0.07MM - CURVED (10)
- #700-431 QWIKSTRIP IPR YELLOW 0.18MM - DOUBLE-SIDED (10)

Kits



- #700-412 QWIKSTRIP IPR RED 0.10MM - SINGLE-SIDED (10)
- #700-422 QWIKSTRIP IPR RED 0.10MM - CURVED (10)
- #700-432 QWIKSTRIP IPR RED 0.20MM - DOUBLE-SIDED (10)

Super Fine



- #700-413 QWIKSTRIP IPR BLUE 0.13MM - SINGLE-SIDED (10)
- #700-423 QWIKSTRIP IPR BLUE 0.13MM - CURVED (10)
- #700-433 QWIKSTRIP IPR BLUE 0.25MM - DOUBLE-SIDED (10)

Fine



- #700-414 QWIKSTRIP IPR GREEN 0.16MM - SINGLE-SIDED (10)
- #700-424 QWIKSTRIP IPR GREEN 0.16MM - CURVED (10)
- #700-434 QWIKSTRIP IPR GREEN 0.30MM - DOUBLE-SIDED (10)

Medium



- #700-415 QWIKSTRIP IPR BLACK 0.18MM - SINGLE-SIDED (10)
- #700-425 QWIKSTRIP IPR BLACK 0.18MM - CURVED (10)
- #700-435 QWIKSTRIP IPR BLACK 0.40MM - DOUBLE-SIDED (10)

Coarse



- #700-436 QWIKSTRIP IPR PURPLE 0.50MM - DOUBLE-SIDED (10)

Extra Coarse



- #700-416 QWIKSTRIP IPR WHITE 0.05MM - SERRATED (10)
- #700-426 QWIKSTRIP IPR WHITE 0.05MM - SERRATED CURVED (10)

XX Coarse



Chain Elastic Dispenser

Provides an efficient way to organize your chain elastics 9 $\frac{7}{8}$ "W x 3 $\frac{3}{8}$ "D x 4 $\frac{1}{4}$ "H.



#700-065

Upright Pliers Rack w/Cover

The compact design holds 12 instruments. The clear cover protects the pliers and easily swings open for easy access. 6" wide x 4" deep x 7" tall.



#700-071

Deluxe Archwire Rack w/Cover

This compact rack has 14 rows of tubes to hold wires securely. It is perfect for right form, global form or D form wires. The clear hinged cover keeps the archwires sanitary and free of dust. 10" wide x 4" deep x 3 $\frac{1}{2}$ " tall



#700-066

Posted/Looped Archwire Rack

An excellent archwire organizer to organize looped and posted wires. Has 12 rows of tubes that are numbered for easy identification. 9 $\frac{1}{2}$ " wide x 4" deep x 2 $\frac{1}{4}$ " tall.



#700-060

Utility Archwire Rack

Features 6 rows of tubes and designed to hold Preformed Utility Archwires. Numbered sizes are: .018 & .016x.016 in 28, 34 and 38mm. 4" wide x 4" deep x 2 $\frac{1}{2}$ " tall



#700-061

Glide-Ties™ Dispenser

An excellent organizer of regular size Glide-Ties™. 18 rows of tubes to hold the most popular ligature colors. The clear hinged cover keeps the ligature ties free of dust particles. Easy access. 10" wide x 6" deep x 3 $\frac{1}{2}$ " tall



#300-226

Glide-Ties™ / Chain Combo Dispenser

Designed to hold 12 each of our most popular colors of matching Glide-Ties and short length chain elastic. The clear hinged cover keeps the elastics free of dust particles. Easy access to pick Glide-Ties™ or cut your chain to length. 7" wide x 9 $\frac{1}{2}$ " deep x 4" tall.



#700-062