




BUCCAL TUBES





- Fluted™ Direct Bond
- Direct Bond
- Weldable
- Specialty Direct Bond
- Light Wire





FLUTED™

DIRECT BOND BUCCAL TUBES

- Color Coded Indent
 - Positioning recesses assure a firm grip during bonding
 - Buccal indent for accurate placement
 - 25% smaller than standard single tubes
 - First molar buccal tube
 - Low profile design and smooth contours for patient comfort
 - Funneled opening helps guide the wire into the buccal tube
 - Compound contoured base to maximize the perfect fit on the buccal surface
 - Second molar buccal tube
- 


BUCCAL TUBES

Maxillary	Description	Torque	Tip	Offset	R/L	.018	.022
	First Molar Single Non-Convertible	-10°	0°	8°	UR ●	#924-210F	#924-310F
		-10°	0°	8°	UL ●	#924-212F	#924-312F
		-14°	0°	8°	UR ●	#924-403F	#924-940F
		-14°	0°	8°	UL ●	#924-404F	#924-942F
	Second Molar Single	-10°	0°	8°	UR ●	#924-033F	#924-133F
		-10°	0°	8°	UL ●	#924-034F	#924-134F
Mandibular	Description	Torque	Tip	Offset	R/L	.018	.022
	First Molar Single Non-Convertible	-20°	0°	5°	LR ●	#924-661F	#924-961F
		-20°	0°	5°	LL ●	#924-662F	#924-962F
		-25°	0°	5°	LR ●	#924-230F	#924-330F
		-25°	0°	5°	LL ●	#924-231F	#924-331F
	Second Molar Single	-20°	0°	5°	LR ●	#924-041F	#924-141F
		-20°	0°	5°	LL ●	#924-042F	#924-142F
		-25°	0°	5°	LR ●	#924-043F	#924-143F
		-25°	0°	5°	LL ●	#924-044F	#924-144F

Maxillary	Description	Torque	Tip	Offset	R/L	.018	.022
	First Molar Single Non-Convertible	0°	0°	0°	UR/LL	#924-201F	#924-301F
		0°	0°	0°	UL/LR	#924-202F	#924-302F
	Second Molar Single	0°	0°	0°	UR/LL	#924-030F	#924-130F
		0°	0°	0°	UL/LR	#924-031F	#924-131F
Mandibular	Description	Torque	Tip	Offset	R/L	.018	.022
	First Molar Single Non-Convertible	0°	0°	0°	LR/UL	#924-202F	#924-302F
		0°	0°	0°	LL/UR	#924-201F	#924-301F
	Second Molar Single	0°	0°	0°	LR/UL	#924-031F	#924-131F
		0°	0°	0°	LL/UR	#924-030F	#924-130F


DIRECT BOND BUCCAL TUBES

BUCCAL TUBES

Maxillary	Description	Torque	Tips	Offset	R/L	.018	.022
	Single Convertible First Molar	0°	0°	0°	UR	#924-801	#924-901
		0°	0°	0°	UL	#924-802	#924-902
		-10°	0°	8°	UR	#924-007	#924-107
		-10°	0°	8°	UL	#924-008	#924-108
	Double Convertible w/.018x.025 Auxiliary Tube	-10°	0°	8°	UR	#924-821F	#924-921
		-10°	0°	8°	UL	#924-822F	#924-922
	Double Convertible w/.045 HG to Gingival	-10°	0°	8°	UR	#924-023G	#924-123G
		-10°	0°	8°	UL	#924-024G	#924-124G
	Double Non-Convertible .045 HG Tube to Gingival	0°	0°	5°	UR		#924-125
		0°	0°	5°	UL		#924-126
	Single Non-Convertible First Molar	0°	0°	0°	UR	#924-301LN	#924-351LN
		0°	0°	0°	UL	#924-302LN	#924-352LN
		-10°	0°	8°	UR	#924-401LN	#924-451LN
		-10°	0°	8°	UL	#924-402LN	#924-452LN
	Solo Second Molar	-10°	0°	0°	UR	#924-405LN	#924-455LN
		-10°	0°	0°	UL	#924-406LN	#924-456LN
Mandibular	Description	Torque	Tips	Offset	R/L	.018	.022
	Single Convertible First Molar	0°	0°	0°	LR	#924-802	#924-902
		0°	0°	0°	LL	#924-801	#924-901
		-20°	0°	0°	LR	#924-863	#924-963
		-20°	0°	0°	LL	#924-864	#924-964
		-25°	0°	5°	LR	#924-021	#924-121
		-25°	0°	5°	LL	#924-022	#924-122
	Double Convertible w/.018x.025 Auxiliary Tube	-20°	0°	5°	LR	#924-663	#924-763
		-20°	0°	5°	LL	#924-664	#924-764
		-25°	0°	5°	LR	#924-523	#924-623
		-25°	0°	5°	LL	#924-524	#924-624
	Double Convertible w/.045 LB to Gingival	-25°	0°	5°	LR	#924-033G	#924-333G
		-25°	0°	5°	LL	#924-034G	#924-334G
	Double Non-Convertible .045 Lip Bumper Tube	0°	0°	5°	LR		#924-335
		0°	0°	5°	LL		#924-336
	Single Non-Convertible First Molar	0°	0°	0°	LR	#924-302LN	#924-352LN
		0°	0°	0°	LL	#924-301LN	#924-351LN
		-20°	0°	0°	LR	#924-206LN	#924-255LN
		-20°	0°	0°	LL	#924-207LN	#924-256LN
		-25°	0°	5°	LR	#924-412LN	#924-462LN
		-25°	0°	5°	LL	#924-413LN	#924-463LN
	Solo Second Molar	-25°	0°	0°	LR	#924-414LN	#924-464LN
		-25°	0°	0°	LL	#924-415LN	#924-465LN

*Individual Buccal Tubes are available in 5/pks.

WELDABLE BUCCAL TUBES

Maxillary	Description	Torque	Tips	Offset	R/L	.018	.022
	Single Convertible First Molar	0°	0°	0°	UR	#942-801	#942-901
		0°	0°	0°	UL	#942-802	#942-902
		-10°	0°	8°	UR	#942-007	#942-107
		-10°	0°	8°	UL	#942-008	#942-108
		-14°	0°	8°	UR	#942-840	#942-940
		-14°	0°	8°	UL	#942-841	#942-941
	Double Convertible w/.018x.025 Auxiliary Tube	-10°	0°	8°	UR	#942-821	#942-921
		-10°	0°	8°	UL	#942-822	#942-922
	Double Convertible w/.045 HG to Gingival	-10°	0°	8°	UR	#942-023G	#942-123G
		-10°	0°	8°	UL	#942-024G	#942-124G
	Double Non-Convertible .045 HG Tube to Gingival	0°	0°	5°	UR		#942-125
		0°	0°	5°	UL		#942-126
	Triple Convertible .045 HG Tube to Occlusal	-10°	0°	8°	UR	#942-823	#942-923
		-10°	0°	8°	UL	#942-824	#942-924
	Single Non-Convertible First and Second Molar	0°	0°	0°	UR	#942-001	#942-101
		0°	0°	0°	UL	#942-002	#942-102
		-10°	0°	8°	UR	#942-211	#942-311
		-10°	0°	8°	UL	#942-212	#942-312
		-14°	0°	8°	UR	#942-842	#942-942
		-14°	0°	8°	UL	#942-843	#942-943
Mandibular	Description	Torque	Tips	Offset	R/L	.018	.022
	Single Convertible	0°	0°	0°	LR	#942-802	#942-902
		0°	0°	0°	LL	#942-801	#942-901
		-20°	0°	5°	LR	#942-861	#942-961
		-20°	0°	5°	LL	#942-862	#942-962
		-25°	0°	5°	LR	#942-021	#942-121
		-25°	0°	5°	LL	#942-022	#942-122
		-12°	0°	0°	LR		#942-931
		-12°	0°	0°	LL		#942-932
	Double Convertible w/.018x.025 Auxiliary Tube	-20°	0°	5°	LR	#942-771	#942-761
		-20°	0°	5°	LL	#942-772	#942-762
		-25°	0°	5°	LR	#942-911	#942-913
		-25°	0°	5°	LL	#942-912	#942-914
	Double Convertible w/.045 Lip Bumper Tube To Gingival	-25°	0°	5°	LR	#942-023G	#942-333G
		-25°	0°	5°	LL	#942-024G	#942-334G
	Double Non-Convertible w/.045 Lip Bumper Tube to Gingival	0°	0°	5°	LR		#942-335
		0°	0°	5°	LL		#942-336
	Single Non-Convertible First and Second Molar	0°	0°	0°	LR	#942-002	#942-102
		0°	0°	0°	LL	#942-001	#942-101
		-20°	0°	5°	LR	#942-213	#942-313
		-20°	0°	5°	LL	#942-214	#942-314
		-25°	0°	5°	LR	#942-231	#942-331
		-25°	0°	5°	LL	#942-232	#942-332


NUSTAR SLB DIRECT BOND TUBE (Passive)



- The easy opening SLB clip converts to an open slot to easily place an archwire
- The clip is then closed to make it a tube
- Low profile
- 80 gauge foil mesh pad

Low profile and round corners enhances patient comfort



Torque in base for precise treatment

80 gauge mesh base provides effective bonding strength





Maxillary	Description	Torque	Offset	R/L	.022
	Direct Bond Molar	-10° -10°	10° 10°	UR UL	#924-511 #924-512
Mandibular	Description	Torque	Offset	R/L	.022
	Direct Bond Molar	-20° -20°	0° 0°	LR LL	#924-513 #924-514

VANISH COPOLYMER CLEAR DIRECT BOND TUBES

Maxillary	Description	Torque	Offset	R/L	.036
	Direct Bond Molar	-10° -10°	10° 10°	UR UL	#924-974 #924-975
Mandibular	Description	Torque	Offset	R/L	.036
	Direct Bond Molar	-20° -20°	0° 0°	LR LL	#924-976 #924-977

BICUSPID DIRECT BOND TUBES

Maxillary	Description	Torque	Tip	.022
	Universal DB Tube - Gingival Hook	-7°	0°	#924-010
Mandibular	Description	Torque	Tip	.022
	Universal DB Tube - Gingival Hook	-17°	0°	#924-012

*Individual Buccal Tubes are available in 5/pks.


LIGHT WIRE BUCCAL TUBES

DIRECT BOND LIGHT WIRE BUCCAL TUBES

Mandibular	Description	Length	Torque	Offset	R/L	.036
	Single w/Hook	4.2mm	0°	0°	UR/LL	#924-000 #924-002
		4.2mm	0°	0°	UL/LR	
	Single w/Hook	6.5mm	0°	0°	UR/LL	#924-000L #924-002L
		6.5mm	0°	0°	UL/LR	

BUCCAL TUBES



WELDABLE LIGHT WIRE BUCCAL TUBES

Mandibular	Description	Length	Torque	Offset	R/L	.036
	Single w/Hook	4.2mm	0°	0°	UR/LL	#924-003 #924-004
		4.2mm	0°	0°	UL/LR	

*Recommended for Conlin and Tip Edge Brackets

*Tip Edge is a registered trademark of TP Orthodontics

DIRECT BOND BUCCAL TUBES WITH POWER HOOKS

Mandibular	Description	Length	Torque	Offset	R/L	.036
	Single w/Hook	4.2mm	-10°	10°	UR	#924-970 #924-971
		4.2mm	-10°	10°	UL	
	Single w/Hook	6.5mm	-20°	0°	LR	#924-972 #924-973
		6.5mm	-20°	0°	LL	

*Individual Buccal Tubes are available in 5/pks.