

# D2 APPLIANCE

## Simplifying Class II Correction

The D2 distalizing appliance was developed to simplify treatment of Class II cases during the initial stage of treatment. The patented D2 design treats to a Class I platform by the intraoral elastic forces. Molars are distal rotated into their desired position. Take advantage of early treatment patient cooperation. Brackets can be placed during or after the Class I correction is completed.



Stainless Steel

## Features & Benefits



One size **fits all** to greatly reduce your inventory



Low profile design for greater **patient comfort**



For **bilateral or unilateral** case treatment

• Simple to use and easy to place

• Patented for both Class II and Class III treatment

• Reduce overall treatment time



**NEW!**  
D2 Appliance Clear

Treat your Class II Aligner cases with fewer Aligners

## 4 EASY STEPS

Etch and prep teeth to be bonded. Bond Posterior Component

Mark & cut wire bar to determined length

Slide anterior component onto wire bar & bond to tooth

Bond lower molar and apply elastic forces.



# D2 APPLIANCE CLEAR

D2 APPLIANCE



Clear



Elastic Protocol

Level 1: ¼" 6½oz  
Level 2: ¼" 8½oz

Available in both Amber and Clear latex-free elastics.

Treat your Class III cases with the same multi-dimensional appliance!



## Elast-O-Loop

A low-profile direct bond molar option to retain elastics.  
#270-209



## Ordering Information

- D2 Appliance Stainless Steel | Set of 2 | #271-600
- D2 Appliance Clear | Set of 2 | #271-600C
- D2 Level 1 Amber Elastics | ¼ in, 6½ oz | #350-311
- D2 Level 2 Amber Elastics | ¼ in, 8½ oz | #350-312
- D2 Level 1 Clear Elastics | ¼ in, 6½ oz | #350-313
- D2 Level 2 Clear Elastics | ¼ in, 8½ oz | #350-314
- 25 Bags of 100 Elastics
- Elast-O-Loop (10/pk) | #270-209