



MOLAR BANDS

- Halo Molar Bands ▀
- Bicuspid Bands ▀
- Plain Bands ▀
- Pre-Welded ▀
- Lingual Attachments ▀

HALO™ MOLAR BANDS



The Halo Molar Bands are anatomically contoured for precise placement and secure fit. The fully tempered molar bands are micro etched inside for maximum retention. Each band is easily identified with a permanent laser ID mark. Available as a plain band or welded with your favorite buccal tube and lingual attachment. Available in sizes 29½ - 44.

MOLAR BANDS

STARTER KIT

Band Size	29½	30	30½	31	31½	32	32½	33	33½	34	34½	35	35½	36	36½	37	37½	38	38½	39	39½	40	40½	41	41½	42	42½	43	43½	44	Total
(Mid) Qty. Right							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1										25
(Mid) Qty. Left							1	1	2	2	2	2	2	2	2	2	2	2	1	1	1										25

INTRODUCTORY KIT

Band Size	29½	30	30½	31	31½	32	32½	33	33½	34	34½	35	35½	36	36½	37	37½	38	38½	39	39½	40	40½	41	41½	42	42½	43	43½	44	Total
(Full) Qty. Right	1	1	1	1	1	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	1	1	1	1	1	1	1	50
(Full) Qty. Left	1	1	1	1	1	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	1	1	1	1	1	1	1	50
(Mid) Qty. Right								4	4	4	4	4	5	5	4	4	4	4	4												50
(Mid) Qty. Left								4	4	4	4	4	5	5	4	4	4	4	4												50

PROFESSIONAL KIT

Band Size	29½	30	30½	31	31½	32	32½	33	33½	34	34½	35	35½	36	36½	37	37½	38	38½	39	39½	40	40½	41	41½	42	42½	43	43½	44	Total
(Full) Qty. Right	1	1	1	1	1	2	2	3	4	4	4	4	5	5	5	4	4	4	4	3	2	2	1	1	1	1	1	1	1	1	75
(Full) Qty. Left	1	1	1	1	1	2	2	3	4	4	4	4	5	5	5	4	4	4	4	3	2	2	1	1	1	1	1	1	1	1	75
(Mid) Qty. Right								6	6	6	6	6	7	7	7	6	6	6	6												75
(Mid) Qty. Left								6	6	6	6	6	7	7	7	6	6	6	6												75

Leatherette Band Box

A great way to store and identify your band inventory

- #900-170UR | #900-170UL
- #900-170LR | #900-170LL



Molar Band Tray w/ Lid

A great way to store and identify your band inventory

- #900-170



BICUSPID BANDS



- Micro etched inside for greater retention.
- Laser etched sizes for easy identification.
- Ideally tempered to adapt easily to the tooth.
- Precisely designed for optimal fit.
- Universal, left, and right sizing.



MAXILLARY BICUSPID KITS INTRODUCTORY KIT

Band Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	Total
(Full) Qty. Right	1	1	1	1	1	2	2	2	3	3	4	4	4	4	3	3	2	2	2	1	1	1	1	1	50
(Full) Qty. Left	1	1	1	1	1	2	2	2	3	3	4	4	4	4	3	3	2	2	2	1	1	1	1	1	50
(Mid) Qty. Right								4	4	4	4	5	5	4	4	4	4	4							50
(Mid) Qty. Left								4	4	4	4	5	5	4	4	4	4	4							50

MANDIBULAR BICUSPID KITS INTRODUCTORY KIT

Band Size	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	Total	
(Full) Qty. Right	1	1	1	1	1	2	2	2	3	3	4	4	4	4	3	3	2	2	2	1	1	1	1	1	1	50
(Full) Qty. Left	1	1	1	1	1	2	2	2	3	3	4	4	4	4	3	3	2	2	2	1	1	1	1	1	1	50
(Mid) Qty. Right								4	4	4	4	5	5	4	4	4	4	4	4							50
(Mid) Qty. Left								4	4	4	4	5	5	4	4	4	4	4	4							50

Weldable Brackets	Section	Torque	Lingual Options	SKU
.022 Upper Bicuspid-Hook	Right	-7 Torque with Hook	Plain, Lug, or Small Cleat	#901-151
	Left			#901-152
.022 Lower Bicuspid-Hook	Right	-17 Torque with Hook	Plain, Lug, or Small Cleat	#901-161
	Left			#901-162









Bicuspid Band Tray

- A great way to store and identify your Bicuspid Band inventory.









[#900-171](#)

WELDABLE LINGUAL ATTACHMENTS

Item	Description	Qty	SKU#
	Lingual Button This low profile button may be welded to your bicuspid or molar band to hold elastics, springs and wires.	30/pk	#270-106
	Lingual Cleat (Small) May be welded to your molar band to help in seating and holding of elastics, springs and wires.	100/pk	#270-214
	Lingual Cleat (Large) May be welded to either a bicuspid or molar band to help with seating and elastic forces. This small molar cleat is welded to all of our molar assemblies unless specified otherwise.	100/pk	#270-213
	Lingual Sheath .036x.072 with latching distal indent.	10/pk	#270-207
	Lingual Sheath with M&D Hook Has a latching Indent	10/pk	#270-207R #270-207L
	Lingual Sheath with Vertical Ball Hook Has a latching Indent	10/pk	#270-201R #210-201L

MOLAR BANDS

DIRECT BOND LINGUAL ATTACHMENTS

Item	Description	Qty	SKU#
	Micro Direct Bond Button Low profile for patient comfort. Anatomical base for better retention. An excellent design to retain elastics and springs.	20/pk	#270-206
	Bite Ramp Stainless steel design and bonded to the lingual of upper anteriors to open a deep bite or advance the mandible.	10/pk	#270-216
	Tongue Trainer The tongue trainer is designed to stop unwanted habits of thumb sucking and tongue thrusting. The 80 gauge bonding base is anatomically contoured to easily bond to the lingual of the upper or lower anterior teeth.	10/pk	#270-205
	MIM Direct Bond Eyelet An excellent appliance to expose impacted teeth.	10/pk	#270-203
	Composite Direct Bond Button Aesthetically pleasing, low profile button, with a dove-tailed, grooved base for a secure bond. Durable, and stain resistant.	10/pk	#270-204
	Ceramic Direct Bond Button Aesthetic button with a low profile, rounded base, and mechanical lock.	10/pk	#270-202
	Eruption Appliance 24K gold plated for nickel sensitive patients, 1" (2.5cm) traction chain 3.5mm diameter contoured base, 80 gauge foil mesh base, with an eyelet welded to base.	5/pk	#270-650
	Elast-O-Loop Direct Bond Hook Low-profile direct bond hook attachment to retain elastics.	10/pk	#270-209