

BRACKET SYSTEMS

Brackets



PINNACLE
STAINLESS STEEL BRACKET SYSTEM

Revolutionary Strength
and Consistency in a
One-Piece Bracket

MADE IN THE
USA

COMPARE TO:
Mini Master™ by American Orthodontics

Intelligent Hook Design
Smooth rounded hooks
reduce soft tissue irritation

Cross Scribe Lines
Unique vertical and horizontal scribe lines in
the center of the bracket head for superior
visual cues. The patented, raised horizontal line
on the bracket base further enhances visibility
for optimal bonding placement accuracy

Micro-Etched Pylon Base Design
Created from the original Marquis™ base,
Pinnacle's redesigned version features angled
pylons that offer increased surface area for
superior retention and bond strength




Permanent Identification System
With a removable color coded dot
and permanently recessed I.D. mark

Rounded Archwire Slot Contours
Mesial and distal slot entry points are
rounded to prevent archwire notching when
rotations are present and during retraction/
protraction procedures

Optimized Compound Contour Base
Brackets are individually contoured for an ideal
fit to each tooth for correct positioning

Identification System



© 2021 Ortho Technology, Inc. Pinnacle is a registered trademark of Ortho Technology. Mini Master is trademark of American Orthodontics.



Brackets

One of the Most Modern Brackets in the Industry

Pinnacle combines classic design with modern elements for better treatment efficiency. Through advanced Metal Injection Molding (MIM) technology, more intricate contours and structures are achieved for reduced friction, optimum patient comfort, and enhanced visual placement marks. A specialized process leaves Pinnacle with a brilliant high-shine surface.

Helps Reduce the Need for Canine Bracket Re-positioning

Doctors struggling with rotations from retraction procedures are able to place Pinnacle more accurately at the beginning of treatment with anti-rotation in the slot.[†]

Pinnacle®

Ortho Technology's version of the Roth Prescription* with Anti-Rotation

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.130"	R	360-91-1100	360-92-1100
				L	360-91-2100	360-92-2100
2 2 Lateral	+8°	+9°	.109"	R	360-91-1200	360-92-1200
				L	360-91-2200	360-92-2200
3 3 Cuspid	-2°	+13°	.117"	R	360-91-1301	360-92-1301
				L	360-91-2301	360-92-2301
3 3 Cuspid with hook	-2°	+13°	.117"	R	360-91-1311	360-92-1311
				L	360-91-2311	360-92-2311
4 5 4 5 Bicuspid	-7°	0°	.120"	R	360-91-1902	360-92-1902
				L	360-91-2902	360-92-2902
4 5 4 5 Bicuspid with hook	-7°	0°	.120"	R	360-91-1912	360-92-1912
				L	360-91-2912	360-92-2912
4 5 4 5 Bicuspid w/offset pad	-7°	0°	.120"	R	360-91-1922	360-92-1922
				L	360-91-2922	360-92-2922
4 5 4 5 Bicuspid w/hook & offset	-7°	0°	.120"	R	360-91-1932	360-92-1932
				L	360-91-2932	360-92-2932

*The Ortho Technology version of the Roth Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

Maxillary Anti-Rotation

Centrals & Laterals = 0°
Cuspids = +4°M
Bicuspids = +2°D

Mandibular Anti-Rotation

Lower Anterior = 0°
Cuspids = +2°M
Bicuspids = +4°D

Low Profile Smooth Facial Contours

For increased patient comfort

Ample Under Tie-Wing Area

Allows for easy and secure ligation

Anti-Rotation in the Slot

[†]Ortho Technology's version of the Roth prescription features anti-rotation in the archwire slot to prevent unwanted rotation during sliding mechanics

Patient Kits

Description

	.018" Slot	.022" Slot
Upper/Lower 5x5	360-91-2000	360-92-2000
Upper/Lower 5x5 Hooks on Cuspids	360-91-2010	360-92-2010
Upper/Lower 5x5 Hooks 3, 4, & 5	360-91-2011	360-92-2011
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-91-2020	360-92-2020
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-91-2030	360-92-2030
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	360-91-2031	360-92-2031

Item Numbers

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 2 1 2 Lower Anterior	-1°	0°	.096"	U	360-91-6800	360-92-6800
3 3 Cuspid	-11°	+7°	.119"	R	360-91-3301	360-92-3301
				L	360-91-4301	360-92-4301
3 3 Cuspid with hook	-11°	+7°	.119"	R	360-91-3311	360-92-3311
				L	360-91-4311	360-92-4311
4 4 1st Bicuspid	-17°	0°	.119"	R	360-91-3402	360-92-3402
				L	360-91-4402	360-92-4402
4 4 1st Bicuspid with hook	-17°	0°	.119"	R	360-91-3412	360-92-3412
				L	360-91-4412	360-92-4412
4 4 1st Bicuspid w/offset pad	-17°	0°	.119"	R	360-91-3422	360-92-3422
				L	360-91-4422	360-92-4422
4 4 1st Bicuspid w/hook & offset	-17°	0°	.119"	R	360-91-3432	360-92-3432
				L	360-91-4432	360-92-4432
5 5 2nd Bicuspid	-22°	0°	.119"	R	360-91-3502	360-92-3502
				L	360-91-4502	360-92-4502
5 5 2nd Bicuspid with hook	-22°	0°	.119"	R	360-91-3512	360-92-3512
				L	360-91-4512	360-92-4512
5 5 2nd Bicuspid w/offset pad	-22°	0°	.119"	R	360-91-3522	360-92-3522
				L	360-91-4522	360-92-4522
5 5 2nd Bicuspid w/hook & offset	-22°	0°	.119"	R	360-91-3532	360-92-3532
				L	360-91-4532	360-92-4532

BRACKET SYSTEMS



Brackets



Pinnacle®

Ortho Technology's version of the McLaughlin, Bennett, & Trevisi Prescription*

Patient Kits Description With Torque on Cuspids

	Item Numbers
Upper/Lower 5x5	360-22-2000
Upper/Lower 5x5 Hooks on Cuspids	360-22-2010
Upper/Lower 5x5 Hooks on 3, 4, 5	360-22-2011
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-22-2020
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-22-2030
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	360-22-2031

Item Numbers .022" Slot

▲ Maxillary Arch	Torq.	Ang.	Width	R/L	.022" Slot
1 1 Central	+17°	+4°	.130"	R	360-22-1100
				L	360-22-2100
2 2 Lateral	+10°	+8°	.109"	R	360-22-1200
				L	360-22-2200
3 3 Cuspid	-7°	+8°	.117"	R	360-22-1300
				L	360-22-2300
3 3 Cuspid with hook	-7°	+8°	.117"	R	360-22-1310
				L	360-22-2310
3 3 Cuspid	0°	+8°	.117"	R	360-22-1340
				L	360-22-2340
3 3 Cuspid with hook	0°	+8°	.117"	R	360-22-1350
				L	360-22-2350
4 5/4 5 Bicuspid	-7°	0°	.120"	R	360-22-1900
				L	360-22-2900
4 5/4 5 Bicuspid with hook	-7°	0°	.120"	R	360-22-1910
				L	360-22-2910
4 5/4 5 Bicuspid with offset	-7°	0°	.120"	R	360-22-1920
				L	360-22-2920
4 5/4 5 Bicuspid w/hook & offset	-7°	0°	.120"	R	360-22-1930
				L	360-22-2930

Single Patient Kit Description

Upper/Lower 5x5

1 kit (20 brackets)

Individual Brackets

1 pack (10 brackets)

Patient Kits Description Without Torque on Cuspids

	Item Numbers
Upper/Lower 5x5	360-22-2040
Upper/Lower 5x5 Hooks on Cuspids	360-22-2050
Upper/Lower 5x5 Hooks on 3, 4, 5	360-22-2051
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	360-22-2060
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	360-22-2061
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	360-22-2062

Item Numbers .022" Slot

▼ Mandibular Arch	Torq.	Ang.	Width	R/L	.022" Slot
1/2 1/2 Lower Anterior	-6°	0°	.096"	U	360-22-6800
3 3 Cuspid	-6°	+3°	.119"	R	360-22-3300
				L	360-22-4300
3 3 Cuspid with hook	-6°	+3°	.119"	R	360-22-3310
				L	360-22-4310
3 3 Cuspid	0°	+3°	.119"	R	360-22-3340
				L	360-22-4340
3 3 Cuspid with hook	0°	+3°	.119"	R	360-22-3350
				L	360-22-4350
4 4 1st Bicuspid	-12°	+2°	.119"	R	360-22-3400
				L	360-22-4400
4 4 1st Bicuspid with hook	-12°	+2°	.119"	R	360-22-3410
				L	360-22-4410
4 4 1st Bicuspid with offset	-12°	+2°	.119"	R	360-22-3420
				L	360-22-4420
4 4 1st Bicuspid w/hook & offset	-12°	+2°	.119"	R	360-22-3430
				L	360-22-4430
5 5 2nd Bicuspid	-17°	+2°	.119"	R	360-22-3500
				L	360-22-4500
5 5 2nd Bicuspid with hook	-17°	+2°	.119"	R	360-22-3510
				L	360-22-4510
5 5 2nd Bicuspid with offset	-17°	+2°	.119"	R	360-22-3520
				L	360-22-4520
5 5 2nd Bicuspid w/hook & offset	-17°	+2°	.119"	R	360-22-3530
				L	360-22-4530

*The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx is not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.



Pinnacle Ideal Occlusion Typodont
Item #: G50-360



WORKS WITH OUR
Votion® Bondable Buccal Tubes

See page 69

